


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/05/16
Time Processed : 11:18:27
Injection Volume : 0.2
Sample Name : Human Papillomavirus (HPV) E7 protein (49 - 57)
Cat No. : RP20249
Lot No. : PE1781610
Theoretical MW : 1120.27
Observed MW : 1119.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH