


Mass Spectrum



Sample Information

Month-Day Processed : 01/27/21
Time Processed : 10:30:10 AM
Injection Volume : 0.6
Sample Name : TAT
Cat No. : RP12784
Lot No. : E0663321A
Theoretical MW : 1558.85
Observed MW : 1558.7

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH