


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 05/23/19
Time Processed : 03:42:20 PM
Injection Volume : 0.3
Sample Name : Angiotensin I, human
Cat No. : RP11251
Lot No. : PE1301905
Theoretical MW : 1296.48
Observed MW : 1296.6

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200