


# Mass Spectrum



Sample Information  
User : Gary  
Month-Day Processed : 02/21/20  
Time Processed : 03:52:55 PM  
Injection Volume : 0.3  
Sample Name : HIV-1 TAT Protein Peptide  
Cat No. : RP10129  
Lot No. : E6525220B  
Theoretical MW : 1559.90  
Observed MW : 1560.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH