


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/22/17
Time Processed : 22:22:19
Injection Volume : 0.2
Sample Name : Cys-TAT(47-57)
Cat. No. : RP20343
Lot. No. : PE4021711
Theoretical MW : 1661.99
Observed MW : 1662.00

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow
T.Flow
B.conc

: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50%H2O/50%MeOH