


Mass Spectrum



Sample Information

Acquired by: Gary
 Month-Day Processed: 05/31/19
 Time Processed: 08:33:38 PM
 Injection Volume: 0.3
 Sample Name: C-Type Natriuretic Peptide (CNP) (1-22), human
 Sample ID: U4419EE040-1
 Theoretical MW: 2197.60
 Observed MW: 2197.5

Interface

Interface: ESI
 Nebulizing Gas Flow: 1.5L/min
 CDL Temp: 250
 Block Temp: 200

Equipment

Equipment: GK11010007
 Interface Bias: +4.5 kV
 Drying Gas Flow: 5 L/min
 T.Flow: 0.2 ml/min
 B.conc: 50% H₂O/50% MeOH