


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/27/20
Time Processed : 9:59:21
Injection Volume : 0.6
Sample Name : TIP 39, Tuberoinfundibular Neuropeptide
Cat No. : RP20322
Lot No. : E5936420C
Theoretical MW : 4504.20
Observed MW : 4504.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH