


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 06/15/18  
Time Processed : 20:04:13  
Injection Volume : 0.2  
Sample Name : Proctolin  
Cat. No. : RP10475  
Lot. No. : PE2501806  
Theoretical MW : 648.75  
Observed MW : 648.45

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias  
Drying Gas Flow  
T.Flow  
B.conc

: +4.5 kV  
: 5 L/min  
: 0.2 ml/min  
: 50% H<sub>2</sub>O/50% MeOH