


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 12/14/18  
Time Processed : 05:40:48 PM  
Injection Volume : 0.1  
Sample Name : Neurotensin  
Cat No. : RP10439  
Lot No. : PE9231812  
Theoretical MW : 1671.92  
Observed MW : 1673.1

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH