


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 10/09/19  
Time Processed : 05:21:58 PM  
Injection Volume : 0.3  
Sample Name : LL-37, acetylated, amidated  
Cat No. : RP20506  
Lot No. : E4942819I  
Theoretical MW : 4534.40  
Observed MW : 4533.9

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

## Equipment

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH