


Mass Spectrum



Sample Information

Month-Day Processed : 11/21/20
Time Processed : 09:12:43 AM
Injection Volume : 0.3
Sample Name : HIF-1 alpha (556-574)
Cat No. : RP20479
Lot No. : E5630620K
Theoretical MW : 2254.60
Observed MW : 2253.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH