


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 09/18/19  
Time Processed : 03:11:33 PM  
Injection Volume : 0.3  
Sample Name : HIV gag peptide (199-207)  
Cat No. : RP20251  
Lot No. : E48503191  
Theoretical MW : 1064.33  
Observed MW : 1064.2

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200