

Sample Name :Lamprey Gonadotropin Releasing Hormone (LGNRH)-III, lamprey
 Cat. No. : RP17322
 Lot. No. : E4900319I
 Time Processed :15:41:09
 Month-Day-Year Processed : 09/23/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
27.00	Pumps	Pump A B.Conc	95	
27.01	Pumps	Pump A B.Conc	5	
35.00	Pumps	Pump A B.Conc	5	
35.00	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0286

Chromatogram



1 Det.A Ch1 / 220nm
 Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	11.095	1559	296	0.023
2	11.240	21573	2734	0.319
3	12.627	14047	2603	0.207
4	12.808	38005	6563	0.561
5	13.020	38253	5365	0.565
6	13.145	12300	2450	0.182
7	13.471	6541573	998312	96.589
8	13.699	105260	21769	1.554
Total		6772570	1040091	100.000