


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/22/19
Time Processed : 07:29:14 PM
Injection Volume : 0.4
Sample Name : Lamprey Gonadotropin Releasing Hormone (LGnRH)-III, lamprey
Cat No. : RP17322
Lot No. : E4900319I
Theoretical MW : 1259.34
Observed MW : 1259.1

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH