


Mass Spectrum



Sample Information

Month-Day Processed : 03/18/21
Time Processed : 9:24:42 AM
Injection Volume : 0.4
Sample Name : Atrial Natriuretic Peptide (ANP) (1-28), human, porcine
Cat No. : RP11927
Lot No. : E1325621C
Theoretical MW : 3080.50
Observed MW : 3080.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH