


Mass Spectrum



Sample Information
Month-Day Processed : 09/08/20
Time Processed : 07:10:49 AM
Injection Volume : 0.4
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No.: RP11110
Lot No.: E38669201
Theoretical MW : 2197.60
Observed MW : 2197.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH