


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/18/19
Time Processed : 02:07:57 PM
Injection Volume : 0.3
Sample Name : HIV p17 Gag (77-85)
Cat No. : RP17316
Lot No. : E4898319I
Theoretical MW : 981.11
Observed MW : 980.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH