


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 03/20/18
Time Processed : 15:57:45
Injection Volume : 0.2
Sample Name : HEL (46-61)
Cat. No. : RP11736 Lot.No.: PE8461803
Theoretical MW : 1753.83
Observed MW : 1754.00

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH