


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/29/18
Time Processed : 01:07:03 PM
Injection Volume : 0.1
Sample Name : AngiotensinII.human
Cat No. : RP11257
Lot No. : PE1381811
Theoretical MW : 1046.18
Observed MW : 1046.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH