


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 12/28/17
 Time Processed : 14:11:05
 Injection Volume : 0.2
 Sample Name : YRGDS Fibronectin Fragment
 Cat. No. : RP20476
 Lot. No. : PE2311712
 Theoretical MW : 596.60
 Observed MW : 596.45

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH