


Mass Spectrum



Sample Information
Month-Day Processed : 11/05/20
Time Processed : 23:24:50
Injection Volume : 0.3
Sample Name : C-telopeptide
Cat No : RP10566
Lot No : E9797120J
Theoretical MW : 868.90
Observed MW : 868.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH