


Mass Spectrum



Sample Information

Month-Day Processed : 11/13/20

Time Processed : 0:00:25

Injection Volume : 0.5

Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human

Cat No. : RP11119

Lot No. : E6026120K

Theoretical MW : 3464.10

Observed MW : 3463.2

Equipment : GK11010007

Interface Bias : +4.5kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2ml/min

B.conc : 50%H2O/50%MeOH

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200