


Mass Spectrum



Sample Information
Month-Day Processed : 10/12/20
Time Processed : 07:06:36 PM
Injection Volume : 0.3
Sample Name : Histatin 5
Cat. No. : RP11240
Lot. No. : E2041720J
Theoretical MW : 3036.29
Observed MW :

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH