


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/10/17
Time Processed : 9:49:58
Injection Volume : 0.2
Sample Name : Neurotensin
Cat No. : RP10439
Lot No. : PE4501708
Theoretical MW : 1671.92
Observed MW : 1672.80

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50% H₂O/50% MeOH