


Mass Spectrum



Sample Information
User : Gary
Month-Day Processed : 01/10/17
Time Processed : 1:00:29
Injection Volume : 1
Sample Name : Enfuvirtide (T-20)
Cat No. : RPI7307
Lot No. : PE1921612
Theoretical MW : 4491.88
Observed MW : 4492.20

Interface
Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias
Drying Gas Flow : +4.5kV
T.Flow : 5 L/min
B.conc : 0.2ml/min
:50%H2O/50%MeOH