


# Mass Spectrum



Sample Information  
Month-Day Processed : 03/26/21  
Time Processed : 9:55:23 PM  
Injection Volume : 0.4  
Sample Name : Jagged-1 (188-204)  
Cat No. : RP20331  
Lot No. : E7133621C  
Theoretical MW : 2107.40  
Observed MW : 2107.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH