


Mass Spectrum



[M+2H]²⁺

[M+3H]³⁺

Sample Information

User : Gary

Month-Day Processed : 07/14/20

Time Processed : 09:43:38 PM

Injection Volume : 0.3

Sample Name : IRBP (1-20), human

Cat No. : RP20269

Lot No. : E9620820G

Theoretical MW : 2194.60

Observed MW : 2194.0

Equipment : GK11010007

Interface Bias : +4.5kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2ml/min

B.conc : 50%H2O/50%MeOH

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200