


Mass Spectrum



Sample Information

Acquired by:	Gary	Equipment	: GK11010007
Month-Day Processed:	07/30/19	Interface Bias	: +4.5 kV
Time Processed:	10:37:12 AM	Drying Gas Flow	: 5 L/min
Injection Volume:	0.3	T.Flow	: 0.2 ml/min
Sample Name:	Ghrelin, rat	B.conc	: 50% H2O/50% MeOH
Cat. No.:	RP10781		
Lot. No.:	E4102419G		
Theoretical MW:	3314.80		
Observed MW:	3315.6		