

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. :RP10769
 Lot No. :E2595121K
 Time Processed :16:06:43
 Month-Day-Year Processed :11/29/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm


Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Chromatogram



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	17.004	31338	5391	1.306
2	17.153	2326052	306223	96.964
3	17.661	24337	2964	1.015
4	17.987	17152	2208	0.715
Total		2398879	316786	100.000