


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 06/14/19
Time Processed : 08:00:28 AM
Injection Volume : 0.3
Sample Name : Angiotensin I, human
Cat. No. : RP11251
Lot. No. : E1187419F
Theoretical MW : 1296.48
Observed MW : 1296.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH