


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 10/19/17
Time Processed : 16:39:37
Injection Volume : 0.2
Sample Name : SLLK, Control Peptide for TSP1 Inhibitor
Sample ID : U5794CJ160-1
Theoretical MW : 459.59
Observed MW : 459.35

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH