


# Mass Spectrum



## Sample Information

Month-Day Processed : 08/05/21  
Time Processed : 8:04:01  
Injection Volume : 0.4  
Sample Name : Calcitonin Gene Related Peptide (CGRP), rat  
Cat No. : RP11095  
Lot No. : E6319921G  
Theoretical MW : 3806.25  
Observed MW : 3806.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH