


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/06/17
Time Processed : 21:46:42
Injection Volume : 1
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : PE5681701
Theoretical MW : 4186.66
Observed MW : 4186.60

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH