


Mass Spectrum



Sample Information

User : Gary  
Month-Day Processed : 10/02/19  
Time Processed : 06:58:08 PM  
Injection Volume : 0.3  
Sample Name : Tetradecapeptide Renin Substrate, human  
Cat No. : RP20216  
Lot No. : E7093519J  
Theoretical MW : 1760.00  
Observed MW : 1759.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH