


# Mass Spectrum



Sample Information  
Month-Day Processed : 11/27/20  
Time Processed : 09:08:31 AM  
Injection Volume : 0.3  
Sample Name : HIV gag peptide (199-207)  
Cat No : RP20251  
Lot No : E5866820K  
Theoretical MW : 1064.33  
Observed MW : 1064.0

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH