


# Mass Spectrum



Sample Information  
Month-Day Processed : 09/03/20  
Time Processed : 06:50:48 PM  
Injection Volume : 0.4  
Sample Name : Cys-TAT(47-57)  
Cat. No. : RP20343  
Lot. No. : E3153120H  
Theoretical MW : 1661.99  
Observed MW : 1661.5

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH