


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/14/20
Time Processed : 05:37:43 PM
Injection Volume : 0.4
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E2352120A
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH