


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/09/19
Time Processed : 03:36:43 PM
Injection Volume : 0.3
Sample Name : TAT 2-4
Cat No. : RP10131
Lot No. : E6973419J
Theoretical MW : 3215.80
Observed MW : 3216.0

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH