


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/03/19
Time Processed : 06:34:43 PM
Injection Volume : 0.4
Sample Name : SLLK, Control Peptide for TSP1 Inhibitor
Cat No. : RP19989
Lot No. : E39172191
Theoretical MW : 459.59
Observed MW : 459.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH