


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/07/17
Time Processed : 17:52:31
Injection Volume : 0.2
Sample Name : Angiotensin II_{human}
Cat No. : RP11257
Lot No. : PE2251712
Theoretical MW : 1046.18
Observed MW : 1046.40

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH