


Mass Spectrum



Sample Information
Month-Day Processed : 08/27/20
Time Processed : 02:33:33 PM
Injection Volume : 0.4
Sample Name : Insulin Receptor (1142-1153)_{pY}(3 sites)
Cat No. : RP20120
Lot No. : E3155420H
Theoretical MW : 1862.68
Observed MW : 1862.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH