


Mass Spectrum



Sample Information

Acquired by: Gary
 Month-Day Processed: 11/16/19
 Time Processed: 02:54:44 PM
 Injection Volume: 0.6
 Sample Name: Biotin-TAT (47-57)
 Cat. No.: RP20408
 Lot. No.: E6827919K
 Theoretical MW: 1786.20
 Observed MW: 1786.5

Interface: ESI
 Nebulizing Gas Flow: 1.5L/min
 CDL Temp: 250
 Block Temp: 200

Equipment: GK11010007
 Interface Bias: +4.5 kV
 Drying Gas Flow: 5 L/min
 T.Flow: 0.2 ml/min
 B.conc: 50% H2O/50% MeOH