


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 03/10/18
 Time Processed : 10:11:44
 Injection Volume : 0.2
 Sample Name : LL-37, scrambled
 Cat. No. : RP20332
 Lot. No. : PE8801803
 Theoretical MW : 4493.30
 Observed MW : 4494.00

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias
 Drying Gas Flow
 T.Flow
 B.conc

: +4.5 kV
 : 5 L/min
 : 0.2 ml/min
 : 50% H₂O/50% MeOH