


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/19/19
Time Processed : 03:33:01 PM
Injection Volume : 0.4
Sample Name : PACAP (1-27), human, ovine, rat
Cat No. : RP10334
Lot No. : E4709019I
Theoretical MW : 3147.62
Observed MW : 3147.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH