


Mass Spectrum



[M+3H]³⁺

[M+2H]²⁺

Sample Information

Month-Day Processed : 09/17/20
Time Processed : 10:16:29 AM
Injection Volume : 0.3
Sample Name : Tetanus toxin (830-843)
Cat No. : RP19980
Lot No. : E7948820H
Theoretical MW : 1611.84
Observed MW : 1611.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH