


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/11/17
Time Processed : 21:00:38
Injection Volume : 0.2
Sample Name : Gastrin-1, human
Cat No. : RPI2740-0.5
Lot No. : PE6151709
Theoretical MW : 2098.21
Observed MW : 2098.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH