

Sample Name :TAT peptide
Cat No. :RP20256
Lot No. :E9327921I
Time Processed :16:58:41
Month-Day-Year Processed :10/07/2021

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
32.00	Pumps	B.Conc	95
32.01	Pumps	B.Conc	5
40.00	Pumps	B.Conc	5
40.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	4.325	8228	1333	0.243
2	5.057	3215873	337665	95.124
3	5.338	105844	8755	3.131
4	5.667	35196	3635	1.041
5	5.955	15585	1330	0.461
Total		3380726	352718	100.000