


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 09/17/19
 Time Processed : 22:18:25
 Injection Volume : 0.3
 Sample Name : Brain Natriuretic Peptide (BNP) (1-32), human
 Cat No. : RP11119
 Lot No. : E4791719I
 Theoretical MW : 3464.10
 Observed MW : 3463.8

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH