


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 06/29/18
Time Processed : 19:20:43
Injection Volume : 0.2
Sample Name : LL-37, scrambled
Cat No. : RP20332
Lot No. : PE2661806
Theoretical MW : 4493.30
Observed MW : 4494.80

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH