

Sample Name :Brain Natriuretic Peptide (BNP) (1-32), rat
 Cat No :RP11121
 Lot No :E9178421C
 Time Processed :15:16:49
 Month-Day-Year Processed :03/27/2021


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	12.004	378332	24927	1.726
2	12.220	20864183	2305062	95.161
3	12.413	682723	131764	3.114
Total		21925238	2461753	100.000